

ANALIZA BIAŁEK



100013526

Laboratory of Analysis

Client: BIO PEPTIDES
Order number: 100013526
Name: BIO PEPTIDES UK CJC-NO DAC/IPA
Number of samples: 1 vial

Batch number: CJ10553
Receive date: 2026-03-20
Start date: 2026-04-01
Enda date: 2026-04-03

Sample identification

Substance	Method of taking a sample	Form
CJC-1295 no DAC, Ipamorelin	Sample provided by customer	White powder in glass vial

Results

Substance	Sample name	Content	Purity
CJC-1295 no DAC	BIO PEPTIDES UK CJC-NO DAC/IPA 5/5MG	5,39mg	99,54%
Ipamorelin	BIO PEPTIDES UK CJC-NO DAC/IPA 5/5MG	4,90mg	99,54%

Research comment

The test sample was applied in the amount of 3 μ l to the BioBasic 8 Thermo Scientific column. During the course of the analysis with the method a linear gradient in the solutions of H₂O and Acetonitrile with the addition of 0.1% trifluoroacetic acid, a peaks of the sought substances appeared. Its concentration, when compared to the reference chromatogram show the table above. Height the measurement error is + - 0,05 mg. The actual state of the sample sent before the test is presented on the photo below

Photo of the test sample



Password

AQ63TR34


Analiza Białek sp. z o.o.
KRS: 0000966168
NIP: 8992922970
REGON: 521751875

Headquarters and Laboratory:
ul. Duńska 9, 54-427 Wrocław
Budynek DELTA

Lab technician:
Analysis and research specialist
mgr Arkadiusz Zajęc



Chromatogram



Method

Time	H2O + 0,1% TFA	ACN + 0,1% TFA
0min	95%	5%
2min	95%	5%
60min	5%	95%
62min	5%	95%

Password

AQ63TR34